

For Research Use Only

PIK3CB Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67121-2-PBS



Basic Information

Catalog Number:

67121-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17505

GenBank Accession Number:

BC114432

GeneID (NCBI):

5291

UNIPROT ID:

P42338

Full Name:

phosphoinositide-3-kinase, catalytic,
beta polypeptide

Calculated MW:

1070 aa, 123 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C6C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

