

For Research Use Only

mFC tag Monoclonal antibody

Catalog Number: 67122-1-Ig



Basic Information

Catalog Number:

67122-1-Ig

Concentration:

1000 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27844

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

26 kDa

Observed MW:

25-30 kDa

Purification Method:

Protein G purification

CloneNo.:

2F5A2

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : Mouse IgG,

Background Information

This antibody detects the mouse Fc tag. It's mouse Fc gamma 2a.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



67122-1-Ig was lysed and subjected to SDS-PAGE followed by western blot with 67122-1-Ig and isotype control antibody 66360-1-Ig at dilution of 1:5000 incubated at room temperature for 1.5 hours. HRP-goat anti-mouse kappa chain secondary antibody was used for detection. .