

For Research Use Only

RGS2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67124-3-PBS



Basic Information

Catalog Number: 67124-3-PBS	GenBank Accession Number: BC007049	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 5997	CloneNo.: 2H9E1
Isotype: IgG1	UNIPROT ID: P41220	
Immunogen Catalog Number: AG28603	Full Name: regulator of G-protein signaling 2, 24kDa	
	Calculated MW: 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

