

For Research Use Only

# EHF Monoclonal antibody

Catalog Number: 67125-1-Ig

Featured Product

1 Publications



## Basic Information

Catalog Number:

67125-1-Ig

Size:

1000  $\mu$ g/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28769

GenBank Accession Number:

BC038995

GeneID (NCBI):

26298

UNIPROT ID:

Q9NZC4

Full Name:

ets homologous factor

Calculated MW:

35 kDa

Observed MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

1C12E4

Recommended Dilutions:

WB 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IHC

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : LNCaP cells, human brain tissue, A431 cells, NCCIT cell, T-47D cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Junqiang Li	33520362	Am J Cancer Res	IHC

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



LNCaP cells were subjected to SDS PAGE followed by western blot with 67125-1-Ig (EHF antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.