

For Research Use Only

RIAM,APBB1IP Monoclonal antibody, PBS Only (Capture)

Catalog Number:67143-2-PBS



Basic Information

Catalog Number: 67143-2-PBS	GenBank Accession Number: BC035636	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 54518	CloneNo.: 2G12G1
Isotype: IgG1	UNIPROT ID: Q7Z5R6	
Immunogen Catalog Number: AG4602	Full Name: amyloid beta (A4) precursor protein- binding, family B, member 1 interacting protein	
	Calculated MW: 666 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51597-1, RIAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67143-2-PBS. Detection antibody: 67143-3-PBS.



Cytometric bead array standard curve of MP51597-1, RIAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67143-2-PBS. Detection antibody: 67143-3-PBS. Standard: Ag4602. Range: 0.781-100 ng/mL