

For Research Use Only

RIAM,APBB1IP Monoclonal antibody, PBS Only (Detector)

Catalog Number:67143-3-PBS



Basic Information

Catalog Number:

67143-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4602

GenBank Accession Number:

BC035636

GeneID (NCBI):

54518

UNIPROT ID:

Q7Z5R6

Full Name:

amyloid beta (A4) precursor protein-
binding, family B, member 1
interacting protein

Calculated MW:

666 aa, 73 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C3F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51597-1, RIAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67143-2-PBS. Detection antibody: 67143-3-PBS.



Cytometric bead array standard curve of MP51597-1, RIAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67143-2-PBS. Detection antibody: 67143-3-PBS. Standard: Ag4602. Range: 0.781-100 ng/mL.