

For Research Use Only

DDX4, VASA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67147-4-PBS



Basic Information

Catalog Number:

67147-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11630

GenBank Accession Number:

BC047455

GeneID (NCBI):

54514

UNIPROT ID:

Q9NQI0

Full Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 4

Calculated MW:

690aa, 76 kDa; 724aa, 79 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E2G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sample test of MP51329-1, DDX4,VASA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67147-3-PBS. Detection antibody: 67147-4-PBS.



Sample test of MP51329-2, DDX4,VASA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67147-1-PBS. Detection antibody: 67147-4-PBS.



Cytometric bead array standard curve of MP51329-1, DDX4,VASA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67147-3-PBS. Detection antibody: 67147-4-PBS. Standard:Ag11630. Range: 0.195-200 ng/mL.



Cytometric bead array standard curve of MP51329-2, DDX4,VASA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67147-1-PBS. Detection antibody: 67147-4-PBS. Standard:null. Range: null.