

For Research Use Only

DHX9 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67153-2-PBS

Basic Information

Catalog Number:

67153-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12104

GenBank Accession Number:

BC014246

GeneID (NCBI):

1660

UNIPROT ID:

Q08211

Full Name:

DEAH (Asp-Glu-Ala-His) box
polypeptide 9

Calculated MW:

1270 aa, 141 kDa

Purification Method:

Protein G purification

CloneNo.:

3B9C9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50161-1, DHX9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67153-2-PBS. Detection antibody: 67153-1-PBS. Standard:Ag12104. Range: 1.563-100 ng/mL.