

For Research Use Only

FABP4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67167-2-PBS



Basic Information

Catalog Number: 67167-2-PBS	GenBank Accession Number: BC003672	Purification Method: Protein G Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 2167	CloneNo.: 2H2D11
Source: Mouse	UNIPROT ID: P15090	
Isotype: IgG1	Full Name: fatty acid binding protein 4, adipocyte	
Immunogen Catalog Number: AG8565	Calculated MW: 132 aa, 15 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51364-1, FABP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67167-2-PBS. Detection antibody: 67167-3-PBS.

Cytometric bead array standard curve of MP51364-1, FABP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67167-2-PBS. Detection antibody: 67167-3-PBS. Standard: Ag8565. Range: 1.563-400 ng/mL