

For Research Use Only

# FABP4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67167-3-PBS



## Basic Information

<b>Catalog Number:</b> 67167-3-PBS	<b>GenBank Accession Number:</b> BC003672	<b>Purification Method:</b> Protein G Magarose purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 2167	<b>CloneNo.:</b> 1A5A7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P15090	
<b>Isotype:</b> IgG1	<b>Full Name:</b> fatty acid binding protein 4, adipocyte	
<b>Immunogen Catalog Number:</b> AG8565	<b>Calculated MW:</b> 132 aa, 15 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP51364-1, FABP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67167-2-PBS. Detection antibody: 67167-3-PBS.

Cytometric bead array standard curve of MP51364-1, FABP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67167-2-PBS. Detection antibody: 67167-3-PBS. Standard: Ag8565. Range: 1.563-400 ng/mL