

For Research Use Only

TRIB2 Monoclonal antibody, PBS Only



Catalog Number: 67171-1-PBS

Basic Information

Catalog Number:

67171-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7501

GenBank Accession Number:

BC002637

GeneID (NCBI):

28951

UNIPROT ID:

Q92519

Full Name:

tribbles homolog 2 (Drosophila)

Calculated MW:

39 kDa

Observed MW:

35-40 kDa

Purification Method:

Protein A purification

CloneNo.:

1A10A6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 67171-1-Ig (TRIB2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67171-1-PBS in a different storage buffer formulation.