

For Research Use Only

ECD Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67179-3-PBS



Basic Information

Catalog Number: 67179-3-PBS	GenBank Accession Number: BC000721	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 11319	CloneNo.: 2D7B6
Isotype: IgG1	UNIPROT ID: O95905	
Immunogen Catalog Number: AG28674	Full Name: ecdysoneless homolog (Drosophila)	
	Calculated MW: 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51594-1, ECD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67179-2-PBS. Detection antibody: 67179-3-PBS.



Cytometric bead array standard curve of MP51594-1, ECD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67179-2-PBS. Detection antibody: 67179-3-PBS. Standard: Ag28674. Range: 1.563-200 ng/mL.