

For Research Use Only

IRAK4 Monoclonal antibody, PBS Only

Catalog Number: 67180-1-PBS



Basic Information

Catalog Number:

67180-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29174

GenBank Accession Number:

BC013316

GeneID (NCBI):

51135

UNIPROT ID:

Q9NWZ3

Full Name:

interleukin-1 receptor-associated
kinase 4

Calculated MW:

52 kDa

Observed MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

1A4C4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67180-1-Ig (IRAK4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67180-1-PBS in a different storage buffer formulation.