

For Research Use Only

JAM-A/Junctional Adhesion Molecule A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67186-2-PBS

Basic Information

Catalog Number:

67186-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25656

GenBank Accession Number:

BC001533

GeneID (NCBI):

50848

UNIPROT ID:

Q9Y624

Full Name:

F11 receptor

Calculated MW:

33 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H3H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50597-1, F11R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67186-1-PBS. Detection antibody: 67186-2-PBS. Standard:Ag25656. Range: 0.195-50 ng/mL