

For Research Use Only

# HSF1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67189-2-PBS



## Basic Information

Catalog Number:

67189-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29061

GenBank Accession Number:

BC014638

GeneID (NCBI):

3297

UNIPROT ID:

Q00613

Full Name:

heat shock transcription factor 1

Calculated MW:

529 aa, 57 kDa

Purification Method:

Protein G purification

CloneNo.:

1A11A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51711-1, HSF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67189-2-PBS. Detection antibody: 67189-3-PBS.



Cytometric bead array standard curve of MP51711-1, HSF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67189-2-PBS. Detection antibody: 67189-3-PBS. Standard: Ag29061. Range: 0.078-10 ng/mL.