

For Research Use Only

KIF20A Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67190-2-PBS



Basic Information

Catalog Number:

67190-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8916

GenBank Accession Number:

BC012999

GeneID (NCBI):

10112

UNIPROT ID:

O95235

Full Name:

kinesin family member 20A

Calculated MW:

890 aa, 100 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E6D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51304-1, KIF20A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67190-2-PBS. Detection antibody: 67190-3-PBS. Standard:Ag8916. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51304-2, KIF20A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67190-4-PBS. Detection antibody: 67190-2-PBS. Standard:Ag8916. Range: 0.195-100 ng/mL