

For Research Use Only

SIN3A Monoclonal antibody, PBS Only



Catalog Number: 67197-1-PBS

Basic Information

Catalog Number:

67197-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28865

GenBank Accession Number:

BC066364

GeneID (NCBI):

25942

UNIPROT ID:

Q96ST3

Full Name:

SIN3 homolog A, transcription regulator (yeast)

Calculated MW:

145 kDa

Observed MW:

145-150 kDa

Purification Method:

Protein G purification

CloneNo.:

2B12E8

Applications

Tested Applications:

WB, IF, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67197-1-Ig (SIN3A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67197-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using SIN3A antibody (67197-1-Ig, Clone: 2B12E8) at dilution of 1:2500 and Multi-rAb CoraLite® Plus 594-Goat Anti-Mouse Recombinant Secondary Antibody (H+L) (Cat.NO. RGAM004). This data was developed using the same antibody clone with 67197-1-PBS in a different storage buffer formulation.