

For Research Use Only

UBE1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67198-2-PBS



Basic Information

Catalog Number: 67198-2-PBS	GenBank Accession Number: BC013041	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 7317	CloneNo.: 1B9E6
Isotype: IgG1	UNIPROT ID: P22314	
Immunogen Catalog Number: AG8920	Full Name: ubiquitin-like modifier activating enzyme 1	
	Calculated MW: 1058 aa, 118 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51632-1, UBE1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67198-2-PBS. Detection antibody: 67198-1-PBS.

Cytometric bead array standard curve of MP51632-1, UBE1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67198-2-PBS. Detection antibody: 67198-1-PBS. Standard: Ag8920. Range: 0.781-100 ng/mL.