

For Research Use Only

EIF4G1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67199-3-PBS

Basic Information

Catalog Number:

67199-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8357

GenBank Accession Number:

BC007788

GeneID (NCBI):

1981

UNIPROT ID:

Q04637

Full Name:

eukaryotic translation initiation factor 4 gamma, 1

Calculated MW:

1600 aa, 176 kDa

Purification Method:

Protein G purification

CloneNo.:

1E7G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50357-1, EIF4G1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67199-2-PBS. Detection antibody: 67199-3-PBS. Standard:Ag8357. Range: 6.25-100 ng/mL