

For Research Use Only

SQLE Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67206-2-PBS



Basic Information

Catalog Number:

67206-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3266

GenBank Accession Number:

BC017033

GeneID (NCBI):

6713

UNIPROT ID:

Q14534

Full Name:

squalene epoxidase

Calculated MW:

574 aa, 64 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3G11E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50872-1, SQLE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67206-2-PBS. Detection antibody: 67206-3-PBS. Standard:Ag3266. Range: 0.195-100 ng/mL.