

For Research Use Only

# WNT10B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67210-3-PBS



## Basic Information

Catalog Number:

67210-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17114

GenBank Accession Number:

BC096353

GeneID (NCBI):

7480

UNIPROT ID:

O00744

Full Name:

wingless-type MMTV integration site family, member 10B

Calculated MW:

389 aa, 43 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G2B6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51624-1, WNT10B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67210-2-PBS. Detection antibody: 67210-3-PBS.



Cytometric bead array standard curve of MP51624-1, WNT10B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67210-2-PBS. Detection antibody: 67210-3-PBS. Standard: Ag17114. Range: 3.12-400 ng/mL.