

For Research Use Only

p130Cas / BCAR1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67215-3-PBS



Basic Information

Catalog Number:

67215-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10604

GenBank Accession Number:

BC062556

GeneID (NCBI):

9564

UNIPROT ID:

P56945

Full Name:

breast cancer anti-estrogen resistance

1

Calculated MW:

870 aa, 93 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B9G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51156-1, p130Cas / BCAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67215-2-PBS. Detection antibody: 67215-3-PBS. Standard:Ag10604. Range: 0.781-12.5 ng/mL



Cytometric bead array standard curve of MP51156-2, p130Cas / BCAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67215-4-PBS. Detection antibody: 67215-3-PBS. Standard:Ag10604. Range: 1.563-100 ng/mL