

For Research Use Only

# GCK Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67216-2-PBS



## Basic Information

<b>Catalog Number:</b> 67216-2-PBS	<b>GenBank Accession Number:</b> BC001890	<b>Purification Method:</b> Protein G purification
<b>Source:</b> Mouse	<b>GeneID (NCBI):</b> 2645	<b>CloneNo.:</b> 1D5G4
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> P35557	
<b>Immunogen Catalog Number:</b> AG8116	<b>Full Name:</b> glucokinase (hexokinase 4)	
	<b>Calculated MW:</b> 52 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51666-1, GSK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67216-2-PBS. Detection antibody: 67216-1-PBS.



Cytometric bead array standard curve of MP51666-1, GSK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67216-2-PBS. Detection antibody: 67216-1-PBS. Standard: Ag8116. Range: 3.125-400 ng/mL.