

For Research Use Only

MMP14 Monoclonal antibody

Catalog Number: 67225-1-Ig



Basic Information

Catalog Number:

67225-1-Ig

Size:

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6062

GenBank Accession Number:

BC064803

GeneID (NCBI):

4323

UNIPROT ID:

P50281

Full Name:

matrix metalloproteinase 14
(membrane-inserted)

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

2F11A11

Applications

Tested Applications:

ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

