

For Research Use Only

MMP14 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67225-2-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
67225-2-PBS	BC064803	Protein G purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	4323	2G6E12
Source:	UNIPROT ID:	
Mouse	P50281	
Isotype:	Full Name:	
IgG1	matrix metalloproteinase 14	
Immunogen Catalog Number:	(membrane-inserted)	
AG6062	Calculated MW:	
	66 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51349-1, MMP14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67225-2-PBS. Detection antibody: 67225-3-PBS.



Cytometric bead array standard curve of MP51349-1, MMP14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67225-2-PBS. Detection antibody: 67225-3-PBS. Standard:Ag6062. Range: 0.098-100 ng/mL.