

For Research Use Only

MMP14 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67225-3-PBS



Basic Information

Catalog Number:

67225-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6062

GenBank Accession Number:

BC064803

GeneID (NCBI):

4323

UNIPROT ID:

P50281

Full Name:

matrix metalloproteinase 14
(membrane-inserted)

Calculated MW:

66 kDa

Purification Method:

Protein G purification

CloneNo.:

2D3B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51349-1, MMP14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67225-2-PBS. Detection antibody: 67225-3-PBS.



Cytometric bead array standard curve of MP51349-1, MMP14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67225-2-PBS. Detection antibody: 67225-3-PBS. Standard:Ag6062. Range: 0.098-100 ng/mL.