

For Research Use Only

TPM4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67244-2-PBS



Basic Information

Catalog Number:

67244-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6947

GenBank Accession Number:

BC037576

GeneID (NCBI):

7171

ENSEMBL Gene ID:

ENSG00000167460

UNIPROT ID:

P67936

Full Name:

tropomyosin 4

Calculated MW:

248 aa, 29 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E11G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51644-1, TPM4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67244-2-PBS. Detection antibody: 67244-1-PBS.

Cytometric bead array standard curve of MP51644-1, TPM4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67244-2-PBS. Detection antibody: 67244-1-PBS. Standard: Ag6947. Range: 1.563-200 ng/mL.