

For Research Use Only

SIRT3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67245-5-PBS



Basic Information

Catalog Number:

67245-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26280

GenBank Accession Number:

BC001042

GeneID (NCBI):

23410

UNIPROT ID:

Q9NTG7

Full Name:

sirtuin (silent mating type
information regulation 2 homolog) 3
(S. cerevisiae)

Calculated MW:

44 kDa

Purification Method:

Protein G purification

CloneNo.:

1D3G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50696-2, SIRT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67245-4-PBS. Detection antibody: 67245-5-PBS. Standard:Ag26280. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50696-3, SIRT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67245-6-PBS. Detection antibody: 67245-5-PBS. Standard:Ag26280. Range: 0.391-100 ng/mL