

For Research Use Only

TMEM230 Monoclonal antibody

Catalog Number: 67247-1-Ig



Basic Information

Catalog Number:

67247-1-Ig

Size:

1500 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG15080

GenBank Accession Number:

BC011990

GeneID (NCBI):

29058

UNIPROT ID:

Q96A57

Full Name:

chromosome 20 open reading frame

30

Calculated MW:

183 aa, 20 kDa

Observed MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

4D3A7

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse, Rat

Positive Controls:

WB : HeLa cells, COLO 320 cells, Caco-2 cells, HT-29 cells, A549 cells, MCF-7 cells, HSC-T6 cells, NIH/3T3 cells, PC-12 cells, HEK-293 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67247-1-Ig (TMEM230 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.