

For Research Use Only

SIRT5 Monoclonal antibody

Catalog Number: 67257-1-Ig **6 Publications**



Basic Information

Catalog Number: 67257-1-Ig	GenBank Accession Number: BC000126	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 23408	CloneNo.: 1F2B2
Isotype: IgG2b	UNIPROT ID: Q9NXA8	Recommended Dilutions: WB: 1:2000-1:10000 IHC: 1:250-1:1000
Immunogen Catalog Number: AG7743	Full Name: sirtuin (silent mating type information regulation 2 homolog) 5 (S. cerevisiae)	
	Calculated MW: 34 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications: WB, IHC, ELISA	Positive Controls:
Cited Applications: WB, IF, IP	WB : HeLa cells, HEK-293 cells, L02 cells, LNCaP cells, K-562 cells
Species Specificity: human	IHC : human liver cancer tissue, human heart tissue
Cited Species: human, bovine	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

SIRT5, a member of sirtuin family that is a mitochondrial matrix NAD(+)-dependent deacetylase and mono-ADP-ribosyltransferases. Sirtuins play important role in various cellular processes, such as aging, gene silencing and metabolism. SIRT5 mainly locates in . It is an efficient protein lysine desuccinylase and demalonylase. During prolonged fasting, it can activate CPS1, the committed and regulated enzyme of the urea cycle, and contribute to the regulation of blood ammonia levels.

Notable Publications

Author	Pubmed ID	Journal	Application
Qichao Huang	35451091	EMBO J	WB
Yiwen Zhang	40481676	Gynecol Endocrinol	WB
Lu Liu	40322223	Genes Dis	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 67257-1-Ig (SIRT5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 67257-1-Ig (SIRT5 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).