

For Research Use Only

# LAMC3-Specific Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67261-2-PBS



## Basic Information

Catalog Number:

67261-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19568

GenBank Accession Number:

NM\_006059

GeneID (NCBI):

10319

UNIPROT ID:

Q9Y6N6

Full Name:

laminin, gamma 3

Calculated MW:

1575 aa, 171 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A4D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51059-1, LAMC3-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67261-1-PBS. Detection antibody: 67261-2-PBS. Standard:Ag19568. Range: 0.195-100 ng/mL