

For Research Use Only

CD5 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67278-2-PBS

Basic Information

Catalog Number: 67278-2-PBS	GenBank Accession Number: BC027901	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 921	CloneNo.: 1B4E5
Source: Mouse	ENSEMBL Gene ID: ENSG00000110448	
Isotype: IgG1	UNIPROT ID: P06127	
Immunogen Catalog Number: AG11097	Full Name: CD5 molecule	
	Calculated MW: 495 aa, 55 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50028-1, CD5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67278-2-PBS. Detection antibody: 67278-3-PBS. Standard: Ag11097. Range: 12.5-100 ng/mL