

For Research Use Only

CD5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67278-3-PBS

Basic Information

Catalog Number:

67278-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11097

GenBank Accession Number:

BC027901

GeneID (NCBI):

921

ENSEMBL Gene ID:

ENSG00000110448

UNIPROT ID:

P06127

Full Name:

CD5 molecule

Calculated MW:

495 aa, 55 kDa

Purification Method:

Protein G purification

CloneNo.:

1A10A5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50028-1, CD5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67278-2-PBS. Detection antibody: 67278-3-PBS. Standard: Ag11097. Range: 12.5-100 ng/mL