

For Research Use Only

KIAA0430 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67281-1-PBS



Basic Information

Catalog Number: 67281-1-PBS	GenBank Accession Number: BC137165	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 9665	CloneNo.: 2D11D7
Source: Mouse	UNIPROT ID: Q9Y4F3	
Isotype: IgG1	Full Name: KIAA0430	
Immunogen Catalog Number: AG18905	Calculated MW: 1742 aa, 193 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51248-1, KIAA0430 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67281-1-PBS. Detection antibody: 67281-2-PBS. Standard: Ag18905. Range: 6.25-200 ng/mL