

For Research Use Only

DLX3 Monoclonal antibody

Catalog Number: 67282-1-Ig



Basic Information

Catalog Number:

67282-1-Ig

Concentration:

1900 ug/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29313

GenBank Accession Number:

BC028970

GeneID (NCBI):

1747

UNIPROT ID:

O60479

Full Name:

distal-less homeobox 3

Calculated MW:

287 aa, 32 kDa

Observed MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

1H8E1

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat, pig

Positive Controls:

WB : pig brain tissue, HEK-293 cells, JAR cells, NCCIT cells, rat brain tissue, mouse brain tissue, mouse cerebellum tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67282-1-Ig (DLX3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.