

For Research Use Only

DLX3 Monoclonal antibody, PBS Only

Catalog Number: 67282-1-PBS



Basic Information

Catalog Number:

67282-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29313

GenBank Accession Number:

BC028970

GeneID (NCBI):

1747

UNIPROT ID:

O60479

Full Name:

distal-less homeobox 3

Calculated MW:

287 aa, 32 kDa

Observed MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

1H8E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67282-1-Ig (DLX3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67282-1-PBS in a different storage buffer formulation.