

For Research Use Only

ALG13 Monoclonal antibody

Catalog Number: 67283-1-Ig



Basic Information

Catalog Number:

67283-1-Ig

Size:

1500 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14955

GenBank Accession Number:

BC005336

GeneID (NCBI):

79868

UNIPROT ID:

Q9NP73

Full Name:

asparagine-linked glycosylation 13 homolog

Calculated MW:

165 aa, 18 kDa

Observed MW:

130-135 kDa

Purification Method:

Protein G purification

CloneNo.:

3B8G3

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : MCF-7 cells, HeLa cells, HEK-293 cells, K-562 cells, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67283-1-Ig (ALG13 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.