

For Research Use Only

# ALG13 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67283-1-PBS



## Basic Information

Catalog Number:

67283-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14955

GenBank Accession Number:

BC005336

GeneID (NCBI):

79868

UNIPROT ID:

Q9NP73

Full Name:

asparagine-linked glycosylation 13 homolog

Calculated MW:

165 aa, 18 kDa

Observed MW:

130-135 kDa

Purification Method:

Protein G purification

CloneNo.:

3B8G3

## Applications

Tested Applications:

WB, Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67283-1-Ig (ALG13 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67283-1-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP50876-1, ALG13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67283-1-PBS. Detection antibody: 67283-2-PBS. Standard: Ag14955. Range: 3.125-200 ng/mL.