

For Research Use Only

# DAAM1 Monoclonal antibody, PBS Only

Catalog Number: 67287-1-PBS



## Basic Information

Catalog Number:

67287-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6809

GenBank Accession Number:

BC064999

GeneID (NCBI):

23002

UNIPROT ID:

Q9Y4D1

Full Name:

dishevelled associated activator of morphogenesis 1

Calculated MW:

123 kDa

Observed MW:

120 kDa

Purification Method:

Protein A purification

CloneNo.:

3F11F11

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Dishevelled-associated activator of morphogenesis 1 (Daam1) is a member of the formin protein family that is implicated in actin assembly. Daam1 also plays a critical role in mediating Wnt/Fz signaling pathway and regulates cell polarity during development and tissue homeostasis. This antibody well recognizes the endogenous DAAM1 protein around 120 kDa.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67287-1-Ig (DAAM1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67287-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded skin cancer slide using 67287-1-Ig (DAAM1 antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67287-1-PBS in a different storage buffer formulation.