

For Research Use Only

# KLF6 Monoclonal antibody, PBS Only

Catalog Number: 67297-1-PBS



## Basic Information

Catalog Number: 67297-1-PBS	GenBank Accession Number: BC000311	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 1316	CloneNo.: 2F8B3
Source: Mouse	UNIPROT ID: Q99612	
Isotype: IgG1	Full Name: Kruppel-like factor 6	
Immunogen Catalog Number: AG6859	Calculated MW: 32 kDa	
	Observed MW: 32-42 kDa	

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
human, mouse

## Background Information

KLF6(krupple-like factor 6) is a zinc finger transcription factor and tumor suppressor with biological activities and transcriptional targets in growing range. It is highly expressed in placenta followed by spleen, thymus, prostate, testis, small intestine and colon. One key target gene of KLF6 in endoglin, which is a homodimeric cell membrane glycoprotein and TGF- $\beta$  auxiliary receptor, owns a pro-angiogenic role in endothelial cells and is also involved in malignant progression. Also KLF6 induces apoptosis in prostate cancer cells through upregulation of ATF3. KLF6 was normal detected as a 46kDa protein, which is larger than translation product(32kDa), because of post-translational modification including both phosphorylation and glycosylation. This antibody raise against the full length of KLF6 gene of human origin, and can recognize both of KLF6, included KLF6(~32kDa), p-KLF6(~35-46kDa). Otherwise, there may be a non-specific band(50kDa) showed in detection

## Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67297-1-Ig (KLF6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67297-1-PBS in a different storage buffer formulation.