

For Research Use Only

# CD107a Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67300-2-PBS



## Basic Information

Catalog Number:

67300-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28869

GenBank Accession Number:

BC021288

GeneID (NCBI):

3916

ENSEMBL Gene ID:

ENSG00000185896

UNIPROT ID:

P11279

Full Name:

lysosomal-associated membrane  
protein 1

Calculated MW:

45 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F11E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50507-1, CD107a Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67300-2-PBS. Detection antibody: 67300-3-PBS. Standard:Ag28869. Range: 6.25-100 ng/mL