

For Research Use Only

Cyclin E Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67302-2-PBS



Basic Information

Catalog Number:	67302-2-PBS	GenBank Accession Number:	BC035498	Purification Method:	Protein G purification
Size:	1 mg/ml	GeneID (NCBI):	898	CloneNo.:	2C2D2
Source:	Mouse	UNIPROT ID:	P24864		
Isotype:	IgG1	Full Name:	cyclin E1		
Immunogen Catalog Number:	AG28418	Calculated MW:	410 aa, 47 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50932-1, Cyclin E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67302-2-PBS. Detection antibody: 67302-3-PBS. Standard: Ag28418. Range: 1.563-100 ng/mL.