

For Research Use Only

ADH4 Monoclonal antibody, PBS Only



Catalog Number: 67303-1-PBS

Basic Information

Catalog Number:

67303-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9884

GenBank Accession Number:

BC022319

GeneID (NCBI):

127

UNIPROT ID:

P08319

Full Name:

alcohol dehydrogenase 4 (class II), pi polypeptide

Calculated MW:

380 aa, 40 kDa

Observed MW:

40-43 kDa

Purification Method:

Protein G purification

CloneNo.:

3D12E9

Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67303-1-Ig (ADH4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67303-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using ADH4 antibody (67303-1-Ig, Clone: 3D12E9) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67303-1-PBS in a different storage buffer formulation.