

For Research Use Only

# FANCI Monoclonal antibody, PBS Only

Catalog Number: 67304-1-PBS



## Basic Information

Catalog Number:

67304-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15153

GenBank Accession Number:

BC004277

GeneID (NCBI):

55215

UNIPROT ID:

Q9NVI1

Full Name:

Fanconi anemia, complementation group I

Calculated MW:

1328 aa, 149 kDa

Observed MW:

149 kDa, 142 kDa

Purification Method:

Protein G purification

CloneNo.:

1B10E6

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67304-1-Ig (FANCI antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67304-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded skin cancer slide using 67304-1-Ig (FANCI antibody) at dilution of 1:2000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67304-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded skin cancer slide using 67304-1-Ig (FANCI antibody) at dilution of 1:2000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67304-1-PBS in a different storage buffer formulation.