

For Research Use Only

# Alpha 2-Antiplasmin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67308-3-PBS



## Basic Information

Catalog Number:

67308-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG4095

GenBank Accession Number:

BC031592

GeneID (NCBI):

5345

UNIPROT ID:

P08697

Full Name:

serpin peptidase inhibitor, clade F  
(alpha-2 antiplasmin, pigment  
epithelium derived factor), member 2

Calculated MW:

491 aa, 55 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A12G5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51218-1, SERPINF2/Alpha 2-Antiplasmin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67308-2-PBS. Detection antibody: 67308-3-PBS. Standard:Ag4095. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51218-2, SERPINF2/Alpha 2-Antiplasmin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67308-4-PBS. Detection antibody: 67308-3-PBS. Standard:Ag4095. Range: 0.781-100 ng/mL