

For Research Use Only

TRIM25 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67314-2-PBS



Basic Information

Catalog Number:

67314-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29325

GenBank Accession Number:

BC016924

GeneID (NCBI):

7706

UNIPROT ID:

Q14258

Full Name:

tripartite motif-containing 25

Calculated MW:

630 aa, 71 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C5F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51519-1, TRIM25 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67314-1-PBS. Detection antibody: 67314-2-PBS.



Cytometric bead array standard curve of MP51519-1, TRIM25 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67314-1-PBS. Detection antibody: 67314-2-PBS. Standard:Ag29325. Range: 1.563-200 ng/mL.