

For Research Use Only

TRAF2 Monoclonal antibody, PBS Only



Catalog Number: 67315-1-PBS

Basic Information

Catalog Number:

67315-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG25387

GenBank Accession Number:

BC032410

GeneID (NCBI):

7186

UNIPROT ID:

Q12933

Full Name:

TNF receptor-associated factor 2

Calculated MW:

56 kDa

Observed MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

2A4D1

Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67315-1-Ig (TRAF2 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67315-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using TRAF2 antibody (67315-1-Ig, Clone: 2A4D1) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67315-1-PBS in a different storage buffer formulation.