

For Research Use Only

ARL2BP Monoclonal antibody

Catalog Number: 67316-1-Ig



Basic Information

Catalog Number:

67316-1-Ig

Size:

1600 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29261

GenBank Accession Number:

BC003087

GeneID (NCBI):

23568

UNIPROT ID:

Q9Y2Y0

Full Name:

ADP-ribosylation factor-like 2 binding protein

Calculated MW:

19 kDa

Observed MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

1G5C2

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : A549 cells, HeLa cells, Jurkat cells, HEK-293 cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.