

For Research Use Only

NICN1 Monoclonal antibody

Catalog Number: 67328-1-Ig



Basic Information

Catalog Number:

67328-1-Ig

Size:

2200 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG6581

GenBank Accession Number:

BC050653

GeneID (NCBI):

84276

UNIPROT ID:

Q9B5H3

Full Name:

nicolin 1

Calculated MW:

24 kDa

Observed MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

5B7G6

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : U-251 cells, NCCIT cells, human testis tissue

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



U-251 cells were subjected to SDS PAGE followed by western blot with 67328-1-Ig (NICN1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.