

For Research Use Only

WDR19 Monoclonal antibody

Catalog Number: 67330-1-Ig



Basic Information

Catalog Number:

67330-1-Ig

Size:

2200 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG4599

GenBank Accession Number:

BC032578

GeneID (NCBI):

57728

UNIPROT ID:

Q8NEZ3

Full Name:

WD repeat domain 19

Calculated MW:

1342 aa, 152 kDa

Observed MW:

145 kDa

Purification Method:

Protein A purification

CloneNo.:

2D1H4

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, NCCIT cells, NIH/3T3 cells, HEK-293 cells, LNCaP cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67330-1-Ig (WDR19 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.